

PRODUCT SPECIFICATION



Cat. No. A-4756_P00
Product Aldol® 515 L-proline amide, hydrochloride, Biosynth Patent: EP 2427431 and US 8940909 and EP1116914

CAS RN **M. F.** C₂₈H₂₉ClN₄O₂
M. W. 489.02 g/mol

<i>Test</i>	<i>Specification</i>
Color	light yellow to yellow to greenish-yellow
Form	powder
Purity (HPLC)	min. 95.0 area-%
TLC (30 mg/mL in methanol; ethyl acetate/pyridine/acetic acid/water 8:2:1:1; UV 254 nm, ninhydrin)	corresponds to reference material
FT-IR Spectrum	corresponds to reference spectrum
Recommended Retest Period:	720 days

Date of issue: June 26, 2016

BIOSYNTH AG

Luca Deriu
